

**BA014**

**WRIGHT'S STAIN**

**FORMULA**

**Reagents:** Prepare freshly working solution before use.

a) Wright stain	1.00 g
Glycerol	50.00 ml
Methanol, absolute	100.00 ml
b) Stock stain solution	4.00ml
Acetone	3.00 ml
Phosphate buffer (1/15M,pH6.5)	2.00 ml
Distilled water	31.00 ml

Both (a) and (b) are mixed in a Coplin jar.

Phosphate buffer (1/15M,pH6.5)

Potassium dihydrogen phosphate, anhydrous	0.663 g
Disodium phosphate, anhydrous	0.256 g
Distilled water	100.00 ml

**Directions:**

1. Air dry blood film.
2. Flood the slide with Wright stain for 3 minutes for fixing.
3. Slowly add buffer of the same quantity as the stain and mix by blowing on the slide.
4. Keep it for 5 minutes.
5. Wash the slide with neutral distilled water and air dry.
6. Observe under oil immersion lens.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** Blood staining

**Quality Control:**

**Appearance:** Blue Purple clear solution

**Microscopy:**

Erythrocytes :	Yellowish red
Polymorpho nuclears:	Dark purple nucleus, reddish lilac granules, pale pink cytoplasm
Eosinophiles :	Blue nuclei, red to orange red granules, blue cytoplasm
Basophiles :	Purple to dark blue nucleus, dark purple granules
Lymphocytes :	Dark purple nuclei , sky blue cytoplasm
Platelates :	Violet to purple granules

**Storage:** Below 30°C.

**Packing:** 200 ml

Refer disclaimer Overleaf

**BIOMARK Laboratories-INDIA**

[www.biomarklabs.com](http://www.biomarklabs.com)

**TECHNICAL SHEET**

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.