BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

BA111	SULPHOSALICYLIC ACID 20%	
Formula		
Reagent:		
	Sulphosalicylic acid20.00 gDistilled water100.00 ml.	
Directions:	 Divide urine sample in two portions. In one portion, add 3 parts of reagent. Allow to stand for 10 minutes and compare with other portion. 	
Precautions :	 For Laboratory Use. Follow proper, established laboratory procedures in handling a disposing of infectious materials. 	and
Use: To estimate protein content from urine by turbidometric method.		
Quality Contro	ol:	
Appearance: Colourless liquid.		
Storage: Below 30°C		
Packing: 100 n	ml bottle	
Disclaimer:		

Page 01 of 01

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.