

BA067

FOUCHET'S REAGENT

Formula

Reagents:

Trichloroacetic acid	25.00 g
Ferric chloride,10% aqueous	10.00 ml
Distilled water	100.00 ml

Directions:

1. Add 1 g of barium chloride to about 10 ml of urine. Mix and filter. Spread out on a filter paper. When it is partly dry, drop a little reagent on the precipitate.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: For testing of bilirubin in urine.

Quality Control:

Appearance: Yellow coloured solution.

Results: Development of green colour indicates the presence of bilirubin in urine.

Storage: Below 30°C.

Packing : 100 ml

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.