

TECHNICAL SHEET

BA068

KOVAC'S INDOLE REAGENT

Formula

Reagents:

p-dimethyl amino benzaldehyde	5.00 g
Amyl alcohol	75.00 ml.
Hydrochloric acid, concentrated	25.00 ml.

Directions:

Add 0.2 to 0.3 ml. of Kovac's reagent to 5 ml. of 18-24 hours old culture of organisms under study. Formation of red colored ring indicates Indole positive test.

Quality Control:

Appearance: Greenish yellow coloured clear solution.

Microbiological testing:

Characteristic reactions observed when reagent is added to growth in Tryptone broth (B1127) incubated for 18 - 24 hours at 35°C - 37°C.

ORGANISM	INDOLE REACTION
Escherichia coli. ATCC 25922	Positive reaction, red ring at the interface of the medium
Enterobacter aerogenes ATCC 13048	Negative reaction, no red ring

Storage: Below 30°C.

Packing: 100 ml

Warning & Precautions :

In Vitro diagnostic Use only. For professional use only. For Laboratory Use. Read the label before opening the container. Wear protective gloves/protective clothing/ eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets

Use: For detection of presence of indole produced by microorganisms due to tryptophan deamination.

Disposal:

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with clinical sample must be decontaminated and disposed of in accordance with current laboratory techniques.

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.