

**BF058 GC SUPPLEMENT W/ANTIBIOTICS**

**FORMULA**

**Ingredients:**

(Per vial, sufficient for 500 ml medium)

Yeast autolysate	5.000 g
Colistin methane sulphonate	3.750 mg
Dextrose	0.750 g
Trimethoprim	2.500 mg
Sodium bicarbonate	0.075 g
Nystatin	6250 Units
Vancomycin	1.500 mg

**Directions:**

Rehydrate the contents of 1 vial aseptically with 15 ml of sterile distilled water. Mix well and add aseptically to 500 ml of sterile, molten, cooled (45 -50°C) GC Agar Base (B1049) / Thayer Martin Medium Base (B333) along With separately prepared Hemoglobin (BF036). Mix well and pour into Sterile petri plates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use: An** antibiotic and enrichment supplement recommended for the selective isolation of pathogenic Neisseria

**Storage:** Between 2 to 8°C.

**Packing :**5 vial

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.