

TECHNICAL SHEET

B1054	HAEMOPHILUS TEST AGAR BASE					
Formula						
Ingredients :		gms/lit.				
Meat infusion from#		300.00				
Casein enzymic hydrolysate		17.50				
Yeast extract		5.00				
Starch		1.50				
Agar		17.00				
#- Equivalent to Beef infusion from						
Final pH (at 25°C) : 7.4 ± 0.2						
Directions :						
Suspend 21.5 grams in 500 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 50°C and aseptically add the rehydrated contents of 1 vial of Haemophilus Growth Supplement (BF064). Mix well and pour into sterile Petri plates.						
Principle :						
Haemophilus influenzae requires complex growth media. Haemophilus Test Agar Base contains meat infusion from and casein acid hydrolysate, which provide essential nutrients to the organisms. Yeast extract serves as a source of B complex vitamins. Starch acts as a protective colloid against toxic substances present in the medium.						
QC Tests – (I) Dehydrated Medium						
Colour :		Cream to yellow				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		7.4 ± 0.2				
Colour (post autoclaving/heating) :		Cream to Light amber				
Clarity (post autoclaving/heating) :		Clear to slightly opalescent				
(III) Q.C. Test Microbiological						
Cultural characteristics observed with added Haemophilus Growth Supplement (BF064) in 5-7% carbon dioxide after an incubation at 35-37°C for 18-24 hours						
MICROORGANISM (ATCC)		GROWTH				
Haemophilus influenzae (49766)		Luxuriant				
Streptococcus pyogenes (19615)		Good – luxuriant				
Neisseria meningitidis (13090)		Good – luxuriant				
Staphylococcus aureus (25923)		Good – luxuriant				
Enterococcus faecalis (29212)		Good – luxuriant				
Precautions :						
1. For Laboratory Use.						
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.						
Limitations :						
1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.						
Use :						
For the susceptibility testing of Haemophilus influenzae.						
Storage :						
Dehydrated medium- below 30°C Prepared medium– Between 2 to 8°C.						
Packing :						
500 gm bottle						
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B1054		43.0 g/l	11.62 L	7.4± 0.2	Haemophilus growth supplement (BF064)	121°C/15min.

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.