

BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

B1187	LACTOBACILLI HETEROFORM SCREEN BROTH(MRS BROTH, MODIFIED)					
Formula						
Ingredients :		gms/lit.				
Dextrose		20.00				
Proteose peptone		10.00				
Yeast extract		5.00				
Sodium acetate		5.00				
2-Phenylethyl alcohol		3.00				
Ammonium citrate		2.00				
Dipotassium phosphate		2.00				
Magnesium sulphate		0.10				
Manganese sulphate		0.05				
Bromocresol green		0.04				
Cycloheximide		0.004				
Final pH (at 25°C) : 4.3± 0.2						
Directions :						
Suspend 47.2 grams in 1000 ml distilled water containing 1 ml polysorbate 80. Mix thoroughly and dispense in tubes containing inverted Durhams tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. If necessary, adjust the pH with glacial acetic acid after sterilization.						
Principle :						
Proteose peptone and dextrose supply nitrogen, carbon and other elements essential for the growth of Lactobacilli. Polysorbate 80 a mixture of oleic esters, supplies fatty acids required by Lactobacilli. Ammonium citrate, sodium acetate, 2-phenylethyl alcohol and cycloheximide inhibit gram-negative organisms, molds and certain gram-positive bacteria. Certain yeasts are also suppressed because of presence of cycloheximide. Bromocresol green is the pH indicator, which under acidic conditions, changes colour from green to yellow.						
QC Tests - (I) Dehydrated Medium						
Colour :		Light yellow to bluish grey				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		4.3 ± 0.2				
Colour (post autoclaving/heating) :		Green				
Clarity (post autoclaving/heating) :		Clear to slightly opalescent solution in tube				
(III) Q.C. Test Microbiological						
Cultural characteristics observed in presence of 5-10% Carbon dioxide(CO ₂) after an incubation at 35-37°C for upto 3 days.						
MICROORGANISM (ATCC)		GROWTH				
Lactobacillus casei (9595)		Luxuriant				
Lactobacillus acidophilus (4356)		Luxuriant				
Lactobacillus fermentum (9338)		Luxuriant				
Lactobacillus plantarum (8014)		Luxuriant				
Precautions :		1. For Laboratory Use.				
		2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
		3. Cycloheximide is very toxic. Avoid skin contact or aerosol formation and inhalation.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		It is recommended for the isolation and cultivation of Lactobacillus species from salad dressings				
Storage :		Dehydrated medium and prepared medium – Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B1187	47.2 g/l		10.593lit	4.3 ± 0.2	polysorbate 80	121°C/15 min

BIOMARK Laboratories-INDIA
www.biomarklabs.com
TECHNICAL SHEET

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.