

BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

B121		BRAIN HEART INFUSION BROTH				
Formula						
Ingredients :		gms/lit.				
Meat brain infusion powder		12.50				
BHI powder		5.00				
Proteose peptone		10.00				
Dextrose (Glucose)		2.00				
Sodium chloride		5.00				
Disodium phosphate		2.50				
Final pH (at 25°C) : 7.4 ± 0.2						
Directions :						
Suspend 37.0 grams in 1000 ml distilled water. Dispense into bottles or tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. For best results, the medium should be used on the day it is prepared, otherwise, it should be boiled or steamed for a few minutes and then cooled before use.						
Principle :						
Proteose peptone, HM infusion powder and BHI powder serve as sources of carbon, nitrogen, essential growth factors, amino acids and vitamins. Dextrose serves as a source of energy. Disodium phosphate helps in maintaining the buffering action of the medium whereas sodium chloride maintains the osmotic equilibrium of the medium. With the addition of 10 % defibrinated sheep blood, it is useful for isolation and cultivation of Histoplasma capsulatum						
QC Tests - (I) Dehydrated Medium						
Colour :		Cream to yellow				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		7.4 ± 0.2				
Colour (post autoclaving/heating) :		Light yellow to medium Amber				
Clarity (post autoclaving/heating) :		Clear				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after 24 -48 hrs at 35-37°C.						
MICROORGANISM (ATCC)		GROWTH				
Neisseria meningitidis (13090)		good-luxuriant				
Streptococcus pneumoniae (6303)		good-luxuriant				
Streptococcus pyogenes (19615)		good-luxuriant				
Staphylococcus aureus (25923)		good-luxuriant				
Candida albicans (10231)		good-luxuriant				
Enterococcus faecalis (29212)		good-luxuriant				
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		For propagation of pathogenic cocci and other fastidious organisms associated with blood culture work and allied pathological investigation.				
Storage :		Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B121		37g/l	13.513L	7.4 ± 0.2	NIL	121°C / 15 minutes

Rev: Feb 2022

BIOMARK Laboratories-INDIA
www.biomarklabs.com
TECHNICAL SHEET

Refer disclaimer Overleaf

Page 01 of 02

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.
The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 02 of 02

Rev: Feb 2022