

TECHNICAL SHEET

B1256	SABOURAUD DEXTROSE BROTH (SABOURAUD LIQUID MEDIUM)					
Formula						
Ingredients : gms/lit.						
Peptone, special 10.00						
Dextrose (Glucose) 20.00						
Final pH (at 25°C) : 5.6 ± 0.2						
Directions :						
Suspend 30.0 grams in 1000 ml distilled water. Heat if necessary, to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.						
Principle :						
Sabouraud Dextrose Broth contains Peptone special or Mycological peptone contains which provide the carbon and nitrogen required for growth of a wide variety of organisms. Dextrose acts as a source of carbohydrates.						
QC Tests - (I)Dehydrated Medium						
Colour :		Cream to light yellow				
Appearance :		Homogeneous Free Flowing powder				
(II)Rehydrated medium						
pH (post autoclaving/heating) :		5.6 ± 0.2				
Colour (post autoclaving/heating) :		Light amber to light yellow				
Clarity (post autoclaving/heating) :		Clear solution in tubes				
(III)Q.C. Test Microbiological						
Cultural characteristics was observed after an incubation at 20-25°C for 3-5 days.						
MICROORGANISM (ATCC)		GROWTH				
Aspergillus niger (16404)		Luxuriant				
Candida albicans (10231)		Luxuriant				
Candida albicans (2091)		Luxuriant				
Saccharomyces cerevisiae (2601)		Good -luxuriant				
Saccharomyces cerevisiae (9763)		Luxuriant				
Escherichia coli (8739)		Luxuriant (inhibited on media with low pH)				
Escherichia coli (25922)		Good -luxuriant				
Escherichia coli (NCTC 9002)		Luxuriant (inhibited on media with low pH)				
Lactobacillus casei (334)		Luxuriant				
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium. 2. Since it is a general purpose medium, bacterial cultures will also grow . 3. Further isolation and biochemical tests should be carried out for confirmation. 4. Avoid overheating a medium with an acidic pH because this often causes a soft medium.				
Use :		For cultivation of yeasts, molds and aciduric organisms.				
Storage :		Dehydrated medium- below 30°C Prepared medium - Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B1256	30g/l	16.666L	5.6 ± 0.2	NIL	121°C / 15 minutes	

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.