BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

B1263	B1263 STREPTOMYCES MEDIUM						
Formula							
Ingredients: gms/lit.							
Glucose			5.00				
L-Glutamic acid	tamic acid						
Monopotassium phosphate			1.00				
Sodium chloride			1.00				
Magnesium sulphate, 7H2O			0.70				
Ferrous sulphate, 7H2O			0.003				
Agar			25.00				
Final pH (at 25°C): 7.0 ± 0.2							
Directions:							
Suspend 36.73 grams of dehydrated medium in 1000 ml. distilled water. Heat to boiling to dissolve the							
medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and							
dispense as desired.							
Principle:							
Glucose is carbon source. L-Glutamic acid is the nitrogen source. The salts provide essential ions for the							
growth of S. kanamyceticus.							
QC Tests - (I)Deh							
Colour:			Off -white to light yellow				
Appearance :			Homogeneous Free Flowing powder				
(II)Rehydrated medium							
PH (post autoclaving/heating) :			7.0 ± 0.2				
			Light amber				
Clarity (post autoclaving/heating):			clear to slightly opalescent gel				
(III)Q.C. Test Microbiological							
	cteristics observe	ed after a	ın incubatior		for upto 7 days.		
MICROORGANISM(ATCC)				GROWTH			
Streptomyces	ces kanamyceticus (12853) luxuriant						
Precautions:	I. For Laboratory Use.						
	2. Follow proper	Follow proper, established laboratory procedures in handling and disposing of					
	infectious materials.						
Limitations : 1. Since the nutritional requirements of organisms vary, some strains						ns may be	
	encountered that fail to grow or grow poorly on this medium.						
Use:	It is recommended for the cultivation and maintenance of Streptomyces kanamyceticus .						
Storage :	Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.						
Packing:	500 gm. bottle						
Product profile:			on	pH (25°C)	Supplement	Sterilization	
		Preparat	ion (500g)				
B1263	36.73 g/l	13	.612 L	7.0 ± 0.2	NIL	121°C / 15 minutes	

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.