BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

| B1264 SUPER BROTH | | | | | | | |
|---|--|---|--|-----------------|-------------------|---------------------|--|
| Formula | | | | | | | |
| Ingredients : gms/lit. | | | | | | | |
| Casein enzymic hydrolysate 35.00 | | | | | | | |
| Yeast extract | | | | 20.00 | | | |
| Sodium chloride | | | 5.00 | | | | |
| Final pH (at 25°C |): Self | | | | | | |
| Directions : | | | | | | | |
| | | | | | | medium completely. | |
| Dispense as desir | ed and sterilize l | oy autocla | aving at 15 ll | bs pressure | (121°C) for 15 mi | nutes. | |
| Principle : | | | | | | | |
| | | | | | | tamin B complex and | |
| other essential growth nutrients. Sodium chloride maintains osmotic equilibrium. Super Broth is | | | | | | | |
| nutritionally rich hence other organisms can also grow in it easily. | | | | | | | |
| QC Tests – (I)Dehydrated Medium | | | | | | | |
| Colour : | | | | Cream to yellow | | | |
| Appearance : | | | Homogeneous Free Flowing powder | | | | |
| (II)Rehydrated medium | | | | | | | |
| PH (post autoclaving/heating) : | | | Self | | | | |
| Colour (post autoclaving/heating) : | | | Light yellow | | | | |
| Clarity (post autoclaving/heating) : | | | clear solution without any precipitate | | | | |
| (III)Q.C. Test Microbiological | | | | | | | |
| Cultural characteristics observed after an incubation at 25-30°C for 18-24 hours. | | | | | | | |
| | | | ROWTH | | | | |
| Escherichia coli (23724) go | | | od-luxuriant | | | | |
| Escherichia coli (25922) go | | | od-luxuriant | | | | |
| Staphylococcu | is aureus (25923 |) go | od-luxuriant | | | | |
| Precautions : | 1. For Laboratory Use. | | | | | | |
| | 2. Follow proper, established laboratory procedures in handling and disposing of | | | | | | |
| | infectious materials. | | | | | | |
| Limitations : | | Since the nutritional requirements of organisms vary, some strains may be | | | | | |
| | encountered that fail to grow or grow poorly on this medium. | | | | | | |
| Use : | It is used for the mass cultivation of Escherichia coli . | | | | | | |
| Storage : | Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C. | | | | | | |
| Packing : | 500 gm. bottle | | | | | | |
| Product profile: | Reconstitution | Quantity Preparat | on ion (500g) | pH (25°C) | Supplement | Sterilization | |
| B1264 | 60.0 g/l | | .33 L | Self | NIL | 121°C / 15 minutes | |
| Disclaimer: | | | | | | | |

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.