

<b>B133</b>	<b>C.L.E.D. AGAR WITH BROMO THYMOL BLUE</b>					
<b>Formula</b>						
<b>Ingredients :</b>			<b>gms/lit.</b>			
Peptic digest of animal tissue			4.00			
Casein enzymic hydrolysate			4.00			
Meat Extract B#			3.00			
Lactose			10.00			
L-cystine			0.128			
Bromo thymol blue			0.02			
Agar			15.00			
#- Equivalent to Beef extract						
Final pH (at 25°C) : 7.3 ± 0.2						
<b>Directions :</b>						
Suspend 36.15 gms in 1000 ml. distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.						
<b>Principle :</b>						
Beef Extract, Peptic digest of animal tissue and Casein enzyme hydrolysate provide the nitrogen, vitamins and amino acids in CLED Agar, L-Cystine is added as a growth supplement for cystine – dependent coliforms. Lactose is included as a carbon source. Organisms capable of fermenting lactose will lower the pH and change the color of the medium from green to yellow. Bromo Thymol Blue is used as a pH indicator.						
<b>QC Tests – (I) Dehydrated Medium</b>						
Colour :			Cream to yellow			
Appearance :			Homogeneous Free Flowing powder			
<b>(II) Rehydrated medium</b>						
pH (post autoclaving/heating) :			7.3 ± 0.2			
Colour (post autoclaving/heating) :			Green			
Clarity (post autoclaving/heating) :			Slightly opalescent			
<b>(III) Q.C. Test Microbiological</b>						
Cultural characteristics observed after 18 –24 hrs at 35-37°C.						
MICROORGANISM (ATCC )		GROWTH		COLOUR OF COLONY		
Escherichia coli (25922)		Luxuriant		Yellow, opaque, center slightly deeper yellow		
Klebsiella pneumoniae (13883 )		Luxuriant		Yellow to whitish blue		
Proteus vulgaris (13315 )		Luxuriant		Blue		
Salmonella typhi ( 6539 )		Luxuriant		Bluish		
Staphylococcus aureus (25923)		Luxuriant		Deep yellow		
Streptococcus faecalis (29212 )		Luxuriant		Slight yellowish or greenish		
<b>Precautions :</b>		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
<b>Limitations :</b>		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
<b>Use :</b>		For isolation, enumeration and identification of urinary pathogens on the basis of lactose fermentation.				
<b>Storage :</b>		Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
<b>Packing :</b>		500 gm bottle				
<b>Product profile:</b>		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
<b>B133</b>		36.15g/l	13.83L	7.3 ± 0.2	NIL	121°C / 15 minutes

Refer disclaimer Overleaf

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.