

B171	DEXTROSE TRYPTONE BROTH					
Formula						
Ingredients :		gms/lit.				
Casein enzymic hydrolysate		10.00				
Dextrose		5.00				
Bromo cresol purple		0.04				
Final pH (at 25°C) : 6.7 ± 0.2						
Directions :						
Suspend 15.04 grams in 1000 ml distilled water. Heat if necessary, to dissolve the medium completely. Dispense and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.						
Principle :						
Casein enzymic hydrolysate provides nutrients to the organisms. Dextrose serves as an energy source while bromo cresol purple is a pH indicator. Acid producing organisms produce yellow coloured medium.						
QC Tests – (I) Dehydrated Medium						
Colour :		Light yellow to greenish yellow				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
PH (post autoclaving/heating) :		6.7 ± 0.2				
Colour (post autoclaving/heating) :		Purple				
Clarity (post autoclaving/heating) :		Clear				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after 48 hrs.at 55°C.						
MICROORGANISM (ATCC)		GROWTH		COLOUR OF MEDIUM		
Bacillus coagulans (8038)		Good to luxuriant		Yellow		
Bacillus stearothermophilus (7953)		Good to luxuriant		Yellow		
Bacillus brevis (8246)		Good to luxuriant*		Yellow		
Key: * = with or without dextrose fermentation.						
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		For enrichment and cultivation of mesophilic and thermophilic organisms in foods.				
Storage :		Dehydrated medium- below 30°C Prepared medium– Between 2 to 8°C.				
Packing :		500 gm bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B171		15.04 g/l	33.244L	6.7 ± 0.2	NIL	121°C / 15 minutes

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.