

TECHNICAL SHEET

B261	MOTILITY SULPHIDE MEDIUM					
Formula						
Ingredients :		gms/lit.				
Proteose peptone		10.00				
Meat Extract B#		3.00				
L-Cystine		0.20				
Ferric ammonium citrate		0.20				
Sodium citrate		2.00				
Sodium chloride		5.00				
Gelatin		80.00				
Agar		4.00				
#- Equivalent to Beef extract						
Final pH (at 25°C) : 7.3 ± 0.2						
Directions :						
Suspend 10.44 grams in 100 ml warm distilled water. Heat to boiling with constant agitation to dissolve the medium completely. Dispense in tubes in 4 ml amounts and sterilize by autoclaving at 115°C (10 lbs pressure) for 15 minutes. Allow the tubed medium to cool in an upright position.						
Principle :						
Proteose peptone, Meat Extract B provides nitrogen compounds, carbon, sulphur and trace elements essential for bacterial growth. L-cystine and ferric ammonium citrate are the H ₂ S indicators. Agar and gelatin preserves an intact stab line.						
QC Tests - (I) Dehydrated Medium						
Colour :		Light yellow				
Appearance :		Homogenous coarse powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		7.3 ± 0.2				
Colour (post autoclaving/heating) :		Yellow				
Clarity (post autoclaving/heating) :		Slightly opalescent				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after 18 -24 hrs.at 35-37°C.						
MICROORGANISM (ATCC)	GROWTH	MOTILITY	H2S	UREASE		
Enterobacter aerogenes (13048)	Luxuriant	+	-	-		
Escherichia coli (8739)	Luxuriant	+	-	-		
Proteus mirabilis (25933)	Luxuriant	±	+	+		
Salmonella typhimurium (14028)	Luxuriant	+	+	-		
Staphylococcus aureus (25923)	Luxuriant	-	-	-		
Shigella sonnei (25931)	Luxuriant	-	-	-		
Key : ± = motility is temperature dependent. It is more pronounced at 20°C and almost absent at 35°C + = positive reaction, motile organism/blackening of medium / reddish purple colouration.						
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		For detection of motility and hydrogen sulphide production by pure cultures.				
Storage :		Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
Packing :		500 gm bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B261		104.4g/l	4.789L	7.3 ± 0.2	NIL	115°C / 15 minutes

Refer disclaimer Overleaf

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.