

**TECHNICAL SHEET**

<b>B283</b>	<b>Phenol red agar base</b>				
<b>Formula</b>					
<b>Ingredients :</b>		<b>gms/lit.</b>			
Proteose peptone		10.00			
Meat Extract B#		1.00			
Sodium chloride		5.00			
Phenol red		0.025			
Agar		15.00			
#- Equivalent to Beef extract					
Final pH (at 25°C) : 7.4 ± 0.2					
<b>Directions :</b>					
Suspend 31.02 grams in 1000 ml purified / distilled water. Add 5-10 grams of carbohydrate as desired. Heat to boiling to dissolve the medium completely. Dispense in tubes or flasks as desired and sterilize by autoclaving at 15 lbs pressure (121° C) for 15 minutes. Allow the tubed media to cool in slanted position to form slants with deep butts. Note: For critical studies, it is recommended to use filter sterilized carbohydrate which is to be incorporated aseptically in sterile medium base					
<b>Principle :</b>					
Proteose Peptone and Meat Extract B provide the carbon and nitrogen required for good growth in a wide variety of organisms. Sodium Chloride maintains the osmotic balance of the medium. Agar is the solidifying agent. Phenol Red serves as a pH indicator, turning from red – orange to yellow when acid is produced during fermentation of the carbohydrate.					
<b>QC Tests – (I) Dehydrated Medium</b>					
Colour :		Light yellow to pink			
Appearance:		Homogeneous Free Flowing powder			
<b>(II) Rehydrated medium</b>					
pH (post autoclaving/heating) :		7.4 ± 0.2			
Colour (post autoclaving/heating) :		Red to orange red			
Clarity (post autoclaving/heating) :		Slightly opalescent			
<b>(III) Q.C. Test Microbiological</b>					
Cultural characteristics observed after 18 - 24 hrs.at 35 -37°C.					
MICROORGANISM (ATCC )	GROWTH	without carbohydrate, (Acid)	without carbohydrate, (Gas)	with dextrose, (Acid)	with dextrose, (Gas)
Alcaligenes faecalis (8750)	Luxuriant	-	-	-	-
Escherichia coli (25922)	Luxuriant	-	-	+	+
Klebsiella pneumoniae (13883 )	Luxuriant	-	-	+	+
Proteus vulgaris (13315)	Luxuriant	-	-	+	+
Salmonella typhimurium (14028)	Luxuriant	-	-	+	+
Shigella flexneri (12022)	Luxuriant	-	-	+	-
*For acid + = Yellow colour					
<b>Precautions :</b>		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.			
<b>Limitations :</b>		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium. 2.Addition of some carbohydrates may result in an acid reaction and hence 0.1N sodium hydroxide can be added dropwise to restore the original colour taking care not to obtain too deep red or cerise colour. 3. When inoculating tubes, stab gently and do not use a loop. Rough stabbing or using a loop to stab may give the false appearance of gas production when mechanical splitting of the medium is what actually occurred.			

Refer disclaimer Overleaf

**BIOMARK Laboratories-INDIA**[www.biomarklabs.com](http://www.biomarklabs.com)**TECHNICAL SHEET**

<b>Use :</b>	As a basal medium to which carbohydrates may be added for use in fermentation studies of microorganisms.				
<b>Storage :</b>	Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
<b>Packing :</b>	500 gm. bottle				
<b>Product profile:</b>	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
<b>B283</b>	31.02g/l	16.12 L	7.4 ± 0.2	Desired carbohydrate	121°C / 15 minutes

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.