

B342	TRYPTOSE BLOOD AGAR BASE					
Formula						
Ingredients :			gms/lit.			
Tryptose			10.00			
Meat Extract B #			3.00			
Sodium chloride			5.00			
Agar			15.00			
# Equivalent to Beef extract						
Final pH (at 25°C) : 7.2 ± 0.2						
Directions :						
Suspend 33 grams in 950 ml purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool the autoclaved medium to 45 - 50°C and aseptically add 5% v/v sterile defibrinated blood. Mix thoroughly, avoiding air bubbles and pour into sterile Petri plates.						
Principle :						
Tryptose is the source of nitrogen, carbon and amino acids in Tryptose Blood Agar Base. Meat extract B provides additional nitrogen, Agar is the solidifying agent. Sodium chloride maintains osmotic balance. Supplementation with 5-10% blood provides additional growth factors for fastidious microorganisms, and is used to determine hemolytic patterns of bacteria.						
QC Tests – (I) Dehydrated Medium						
Colour :		Cream to yellow				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		7.2 ± 0.2				
Colour (post autoclaving/heating) :		a) Basal medium : Yellow b) With addition of 5% v/v defibrinated sterile blood : Cherry red.				
Clarity (post autoclaving/heating) :		a) Clear to slightly opalescent b) Opaque				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after 18 – 48 hrs. at 35-37°C.						
	MICROORGANISM (ATCC)	GROWTH W/O BLOOD	GROWTH W BLOOD	Haemolysis		
	Neisseria meningitidis (13090)	Good - luxuriant	Luxuriant	None		
	Staphylococcus aureus (25923)	Good - luxuriant	Luxuriant	Beta/Gamma		
	Staphylococcus epidermidis (12228)	Good - luxuriant	Luxuriant	Gamma		
	Streptococcus pneumoniae (6303)	Fair-Good	Good	Alpha		
	Streptococcus pyogenes (19615)	Fair-Good	Good	Beta		
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use:		Recommended for the isolation of fastidious organisms and determining the haemolytic reactions.				
Storage :		Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B342		33g/l	15.15L	7.2 ± 0.2	5% v/v sterile defibrinated blood.	121°C /15 min.

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.