

<b>B374</b>	<b>S.F.P. AGAR BASE</b>					
<b>Formula</b>						
<b>Ingredients :</b>		<b>gms/lit.</b>				
Tryptose		15.00				
Papaic digest of soyabean meal		5.00				
Yeast extract		5.00				
Sodium bisulphite		1.00				
Ferric ammonium citrate		1.00				
Agar		20.00				
Final pH (at 25°C) :		7.6± 0.2				
<b>Directions :</b>						
Suspend 23.5 gms.in 475ml. distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 50°C.Add 25 ml. of Egg Yolk Emulsion (BF003) and reconstituted contents of 1 vial of S.F.P. Supplement (BF090).Mix well before pouring into sterile petri plates.						
<b>Principle :</b>						
SFP Agar Base contains Tryptose and Papaic digest of soyabean meal as sources of carbon, nitrogen, vitamins and minerals. Yeast extract supplies B –complex vitamins which stimulate bacterial growth. Ferric Ammonium Citrate and Sodium bisulfite are H <sub>2</sub> S indicators. Clostridia reduce sulfite to sulfide, which reacts with iron to form a black iron sulfide precipitate. Polymyxin B and Kanamycin; are inhibitors to organisms other than Clostridium spp. Egg yolk enrichment 50% provides egg yolk lecithin which some clostridia hydrolyze. Agar is the solidifying agent.						
<b>QC Tests – (I)Dehydrated Medium</b>						
	Colour :	Beige to yellow				
	Appearance :	Homogeneous Free Flowing powder				
<b>(II)Rehydrated medium</b>						
	pH (post autoclaving/heating) :	7.6 ± 0.2				
	Colour (post autoclaving/heating) :	a) Basal medium : Amber b) After addition of Egg Yolk Emulsion : Pale yellow				
	Clarity (post autoclaving/heating) :	a) Slightly opalescent b) Opaque				
<b>(III)Q.C. Test Microbiological</b>						
Cultural characteristics observed after 40 – 48 hrs.at 35- 37°C, when incubated anaerobically.						
	MICROORGANISM (ATCC )	GROWTH	COLOUR OF COLONY	LECITHINASE		
	Clostridium perfringens (12924)	Luxuriant	Black	+		
	Escherichia coli (25922)	Inhibited	-	-		
	Key: + = halo or clearing around colony.					
<b>Precautions :</b>		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
<b>Limitations :</b>		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
<b>Use :</b>		For the presumptive identification and enumeration of Clostridium perfringens in foods.				
<b>Storage :</b>		Dehydrated medium- below 30°C Prepared medium– Between 2 to 8°C.				
<b>Packing :</b>		500 gm. bottle				
<b>Product profile:</b>		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
<b>B374</b>		47g/l	10.638L	7.6 ± 0.2	Egg Yolk Emulsion (BF003) and S.F.P. Supplement (BF090)	121°C / 15 minutes

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.