

<b>B464</b>	<b>DEXTROSE TRYPTONE AGAR</b>				
<b>Formula</b>	<b>gms/lit.</b>				
<b>Ingredients :</b>					
Casein enzymic hydrolysate	10.00				
Dextrose	5.00				
Bromo cresol purple	0.04				
Agar	15.00				
Final pH (at 25°C) : 6.7 ± 0.2					
<b>Directions :</b>					
Suspend 30.0 grams in 1000 ml distilled water. Heat to boiling the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.					
<b>Principle :</b>					
Dextrose Tryptone Agar contains Casein enzymic hydrolysate to provide carbon and nitrogen sources for general growth requirements. Dextrose is the carbohydrate source. Bromo Cresol Purple is the pH indicator. Acid producing organisms produce yellow colonies. Agar is the solidifying agent.					
<b>QC Tests – (I) Dehydrated Medium</b>					
Colour :	Light yellow to greenish yellow				
Appearance :	Homogeneous Free Flowing powder				
<b>(II) Rehydrated medium</b>					
PH (post autoclaving/heating) :	6.7 ± 0.2				
Colour (post autoclaving/heating) :	Purple				
Clarity (post autoclaving/heating) :	Clear to slightly opalescent				
<b>(III) Q.C. Test Microbiological</b>					
Cultural characteristics observed after an incubation at 54-56°C for 36-48 hours					
MICROORGANISM (ATCC )	GROWTH		COLOUR OF MEDIUM		
Bacillus coagulans (8038)	Good to luxuriant		Yellow		
Bacillus stearothermophilus (7953)	Good to luxuriant		Yellow		
Bacillus brevis (8246)	Good to luxuriant*		Yellow		
Key: * = with or without dextrose fermentation.					
<b>Precautions :</b>	1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
<b>Limitations :</b>	1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
<b>Use :</b>	For detection and enumeration of mesophilic and thermophilic aerobic organisms in foods.				
<b>Storage :</b>	Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
<b>Packing :</b>	500 gm bottle				
<b>Product profile:</b>	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
<b>B464</b>	30.0 g/l	16.644 L	6.7 ± 0.2	NIL	121°C / 15 minutes

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.