## **BIOMARK Laboratories-INDIA**

## www.biomarklabs.com

## **TECHNICAL SHEET**

B541 GL	UCOSE CITRATE BROTH BASE					
Formula						
Ingredients:		gms,				
Meat Extract B#		5.0	0			
Peptic digest of a	nimal tissue	10.0				
Sodium chloride		5.0	-			
Glucose	1.00					
Sodium citrate 2.00						
# Equivalent to Beef extract						
Final pH (at 25°C): 7.6 <u>+</u> 0.2						
Directions:						
Suspend 23 grams in 1000 ml distilled water. Heat if necessary to dissolve the medium completely.						
Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Dispense as desired.						
Principle:						
Peptic digest of animal tissue and meat extract B provides nitrogenous compounds and other essential						
growth nutrients. Glucose provides the energy sources and sodium citrate as the carbon source.						
	nydrated Medium					
Colour:			Yellow			
Appearance :			Homogeneous Free Flowing powder			
(II)Rehydrated medium						
			$7.6 \pm 0.2$			
Colour (post autoclaving/heating):			Light yellow			
			Clear			
(III)Q.C. Test Microbiological						
Cultural characteristics observed after 18 –24 hrs. at 35-37°C.						
MICROORGANISM (ATCC )			GROV	VTH		
Enterobacter aerogenes (13048)			Luxur	iant		
Escherichia coli (25922)			Luxur	iant		
Precautions: 1. For Laboratory Use.						
2. Follow proper, established laboratory procedures in handling and disposing						ng and disposing of
infectious materials.						
Limitations:  1. Since the nutritional requirements of organisms vary, some strains made encountered that fail to grow or grow poorly on this medium.						ome strains may be
						•
Use :	It is used for cultivation of fastidious microorganisms.					
Storage :	Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.					
Packing:	500 gm bottle					
Product profile:	Reconstitution Quantity on		1	pH (25°C)	Supplement	Sterilization
_		Preparation				
B541	23.0g/l	21.7	'39 L	$7.6 \pm 0.2$	nil	118°C / 15 minutes
	_					

## Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.