

TECHNICAL SHEET

| | | | |
|--|--|---|----------------------|
| B544 | GLUCOSE STARCH AGAR | | |
| Formula | | | |
| Ingredients : | | gms/lit. | |
| Proteose peptone | | 15.00 | |
| Dextrose | | 10.00 | |
| Starch, soluble | | 5.00 | |
| Sodium chloride | | 5.00 | |
| Disodium hydrogen phosphate | | 3.00 | |
| Gelatin | | 20.00 | |
| Agar | | 10.00 | |
| Final pH (at 25°C) : 7.2 ± 0.2 | | | |
| Directions : | | | |
| Suspend 68.0 grams in warm 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Dispense in tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 30 minutes. Allow the tubed medium to cool in an upright position. | | | |
| Principle : | | | |
| Clostridial species are one of the major causes of food poisoning/ gastro-intestinal illnesses. They are gram-positive, spore-forming rods that occur naturally in the soil. Among the family are: <i>Clostridium botulinum</i> , which produces one of the most potent toxins in existence; <i>Clostridium tetani</i> , causative agent of tetanus; and <i>Clostridium perfringens</i> , commonly found in wound infections and diarrhoea cases. The use of toxins to damage the host is a method deployed by many bacterial pathogens. The major virulence factor of <i>C.perfringens</i> is the CPE enterotoxin, which is secreted upon invasion of the host gut, and contributes to food poisoning and other gastrointestinal illnesses. Glucose Starch Agar is used as a basal medium, which with the addition of raffinose, salicin and phenol red indicator is used for detecting <i>C.perfringens</i> . This medium is also recommended by APHA. | | | |
| The medium contains proteose peptone, which supplies the nitrogenous nutrients for <i>C.perfringens</i> . Dextrose is the fermentable carbohydrate source and is fermented by most Clostridia. However, raffinose and salicin are fermented with acid and gas production by only some strains of <i>C.perfringens</i> . Dispense the medium in different tubes and add a few drops of phenol red, the pH indicator, which turns yellow at acidic pH. Gas production is indicated by bubble formation. Gelatin is liquefied by <i>C. perfringens</i> within 48 hours. Sodium chloride maintains the osmotic balance of the medium. | | | |
| QC Tests - (I)Dehydrated Medium | | | |
| Colour : | | Cream to Beige | |
| Appearance : | | Homogeneous Free Flowing powder | |
| (II)Rehydrated medium | | | |
| pH (post autoclaving/heating) : | | 7.2 ± 0.2 | |
| Colour (post autoclaving/heating) : | | Light amber | |
| Clarity (post autoclaving/heating) : | | Clear to slightly opalescentgel forms in tubes as butts | |
| (III)Q.C. Test Microbiological | | | |
| Cultural characteristics observed after an incubation at 35-37°C for 24-72 hours. Dextrose fermentation is detected using phenol red indicator. | | | |
| MICROORGANISM (ATCC) | GROWTH | SALICIN (24 HOURS) | RAFFINOSE (72 HOURS) |
| <i>Clostridium perfringens</i> (12924) | Luxuriant | - | A |
| <i>Clostridium paraperfringens</i> | Luxuriant | AG | - |
| <i>Escherichia coli</i> (25922) | Luxuriant | - | - |
| Key : A = Acid production AG = Acid and Gas production | | | |
| Precautions : | 1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials. | | |

BIOMARK Laboratories-INDIAwww.biomarklabs.com**TECHNICAL SHEET**

| | | | | | |
|-------------------------|---|--------------------------------|-----------|------------|--------------------|
| Limitations : | 1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium. | | | | |
| Use : | Used as a basal medium with the addition of salicin,raffinose and phenol red for detection of Clostridium perfringens. | | | | |
| Storage : | Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C. | | | | |
| Packing : | 500 gm bottle | | | | |
| Product profile: | Reconstitution | Quantity on Preparation (500g) | pH (25°C) | Supplement | Sterilization |
| B544 | 68g/l | 7.35 L | 7.2 ± 0.2 | nil | 121°C / 15 minutes |

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

REV: Jan. 2023