

BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

B578	MIDDLEBROOK 7H11 AGAR BASE WITHOUT MALACHITE GREEN	
Formula		
Ingredients :		gms/lit.
Casein enzymic hydrolysate		1.00
Ammoniumsulphate		0.50
Monopotassium Phosphate		1.50
Disodium phosphate		1.50
Sodium citrate		0.40
Magnesium sulphate		0.05
L-Glutamic acid		0.50
Ferric ammonium citrate		0.04
Pyridoxine		0.001
Biotin		0.0005
Agar		15.00
Final pH (at 25°C) : 6.6± 0.2		
Directions :		
Suspend 10.25 gms. in 450ml distilled water containing 2.5ml glycerol. Boil to dissolve the medium completely. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes. Cool to 50°C. Aseptically add 1 vial of Middlebrook OADC Growth Supplement. Mix thoroughly before dispensing. Add Malachite green if desired.		
Principle :		
The media consists of many inorganic salts which help for growing Mycobacteria. Citric acid formed from sodium citrate helps in retaining inorganic cations in solution. Glycerol supplies carbon and energy. OADC Supplement contains oleic acid, bovine albumin, sodium chloride, dextrose and catalase. Oleic acid and other long chain fatty acids are metabolized by Mycobacteria. Dextrose is an energy source. Catalase neutralizes toxic peroxides, while albumin protects tubercle bacilli from toxic agents. Malachite green partially inhibits other bacteria.		
QC Tests - (I) Dehydrated Medium		
Colour :		Yellow
Appearance :		Homogeneous Free Flowing powder
(II) Rehydrated medium		
pH (post autoclaving/heating) :		6.6 ± 0.2
Colour (post autoclaving/heating) :		Very light amber
Clarity (post autoclaving/heating) :		Slightly opalescent
(III) Q.C. Test Microbiological		
Cultural characteristics observed after 2-4 weeks at 35-37°C.		
MICROORGANISM (ATCC)		GROWTH
Mycobacterium tuberculosis H37 RV (25618)		Good - Luxuriant
Mycobacterium smegmatis (14468)		Good - Luxuriant
Mycobacterium fortuitum (6841)		Good - Luxuriant
Precautions :	1. For Laboratory Use.	
	2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.	

Refer disclaimer Overleaf

BIOMARK Laboratories-INDIAwww.biomarklabs.com**TECHNICAL SHEET**

Limitations :	1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
	2. Negative culture results do not rule out active infection by mycobacteria. Some factors that are responsible for unsuccessful cultures are ; <ul style="list-style-type: none">• The specimen was not representative of the infectious material, i.e. saliva instead of sputum.• The mycobacteria were destroyed during digestion and decontamination of the specimen.• Gross contamination interfered with the growth of the mycobacteria.• Proper aerobic conditions and increased CO₂ tension were not provided during incubation.				
Use :	For isolation, cultivation and sensitivity testing of Mycobacteria.				
Storage :	Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
Packing :	500 gm. bottle				
Product profile:	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B578	20.49g/l	24.40L	6.6 ± 0.2	Middlebrook OADC Growth Supplement Malachite green	121 ⁰ C / 15 minutes

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.