

BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

B623	MINIMAL AGAR DAVIS					
Formula						
Ingredients :		gms/lit.				
Dextrose		1.00				
Dipotassium phosphate		7.00				
Monopotassium phosphate		2.00				
Sodium citrate		0.50				
Magnesium sulphate		0.10				
Ammonium sulphate		1.00				
Agar		15.00				
Final pH (at 25°C) : 7.0 ± 0.2						
Directions :						
Suspend 26.6 gms. in 1000 ml. distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.						
Principles :						
Minimal Agar Davis contain citrate and phosphates as buffers. Ammonium Sulfate is the carbon source. Magnesium is a cofactor for many metabolic reactions. Minimal Agar Davis contains Dextrose as the carbohydrate energy source. Agar is the solidifying agent.						
QC Tests – (I) Dehydrated Medium						
Colour :		Cream to yellowish white				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		7.0 ± 0.2				
Colour (post autoclaving/heating) :		Medium amber				
Clarity (post autoclaving/heating) :		Slightly opalescent				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after 18- 24 hours at 35-37°C.						
MICROORGANISM (ATCC)		GROWTH				
Escherichia coli (13762)		Luxuriant				
Escherichia coli (23724)		Luxuriant				
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		For the isolation and characterization of nutritional mutants of E.Coli.				
Storage :		Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B623	26.6g/l	18.79L	7.0 ± 0.2	NIL	121°C / 15 minutes	

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Rev: December 2020