BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

B624	MINIMAL BRO	MINIMAL BROTH DAVIS					
Formula							
		gm	gms/lit.				
Dextrose		1.00	1.00				
Dipotassium phosphate		7.00	7.00				
Monopotassium phosphate		2.00	2.00				
Sodium citrate		0.50	0.50				
Magnesium sulphate		0.10	0.10				
Ammonium sulphate		1.00	1.00				
Final pH (at 25	ю́С): 7.0 <u>+</u> 0.2						
Directions :							
Suspend 11.6 gms. in 1000 ml. distilled water.Boil to dissolve the medium completely. Sterilize by							
autoclaving at 15 lbs pressure (121°C) for 15 minutes.							
Principles :							
Minimal Broth Davis contain citrate and phosphates as buffers. Ammonium Sulfate is the carbon source.							
Magnesium is a cofactor for many metabolic reactions. Minimal Agar Davis contains Dextrose as the							
carbohydrate energy source. Agar is the solidifying agent.							
QC Tests – (I)Dehydrated Medium							
Colour :			Cream to yellowish white				
Appearance :			Homogeneous Free Flowing powder				
(II)Rehydrated medium							
pH (post autoclaving/heating) :			7.0 ± 0.2				
Colour (post autoclaving/heating) :			Medium amber				
Clarity (post autoclaving/heating) :			Slightly opalescent				
(III) Q.C. Test Microbiological							
Cultural cha	aracteristics observ	ed after 1	8- 24 hours	at 35-37°C.			
MICROORGANISM (ATCC)			ROWTH				
Escherichia coli (13762)		Lu	ixuriant				
Escherichia coli (23724)		Lu	ixuriant				
Precautions :	1. For Laborate	1. For Laboratory Use.					
		2. Follow proper, established laboratory procedures in handling and disposing of					
		infectious materials.					
Limitations :	1. Since the nu	1. Since the nutritional requirements of organisms vary, some strains may be					
		encountered that fail to grow or grow poorly on this medium.					
Use :	For the isolation and characterization of nutritional mutants of E.Coli.						
Storage :							
Packing : 500 gm. bottle							
	e: Reconstitution Quantity Preparat		on ion (500g)	pH (25°C)	Supplement	Sterilization	
B624	11.6g/l		3.10L	7.0 ± 0.2	NIL	121ºC / 15 minutes	

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.