

B626	MODIFIED ROGOSA AGAR (M16 AGAR)					
Formula						
Ingredients :		gms/lit.				
Papaic digest of soyabean meal		5.00				
Tryptose		5.00				
Meat extract B#		5.00				
Yeast extract		2.50				
Dextrose		5.00				
Ascorbic acid		0.50				
Sodium acetate		3.00				
Agar		10.00				
#- Equivalent to Beef extract						
Final pH (at 25°C) : 7.2± 0.2						
Directions :						
Suspend 36 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.						
Principle :						
Papaic digest of soyabean meal, tryptose and Meat extract B provide the essential nutrients like amino acids, minerals etc. Yeast extract supplies vitamin B complex to the lactic streptococci. Dextrose is the fermentable carbohydrate and energy source. Sodium acetate inhibits other contaminating bacteria and suppresses swarming growth. Ascorbic acid provides vitamin C to the organisms.						
QC Tests – (I) Dehydrated Medium						
Colour :		Cream to yellow				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		7.2 ± 0.2				
Colour (post autoclaving/heating) :		Light amber				
Clarity (post autoclaving/heating) :		Clear to slightly opalescent gel				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours in CO ₂ enriched atmosphere						
MICROORGANISM (ATCC)		GROWTH				
Lactobacillus lactis (8000)		good-luxuriant				
Streptococcus cremoris(19257)		good-luxuriant				
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		It is recommended for cultivation and enumeration of lactic streptococci used in manufacture of cheddar cheese.				
Storage :		Dehydrated medium and prepared medium– Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B626		36.0 g/l	13.888 L	7.2 ± 0.2	NIL	121°C / 15 minutes

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET