

**TECHNICAL SHEET**

<b>B704</b>	<b>SA AGAR BASE</b>				
<b>Formula</b>					
<b>Ingredients :</b>		<b>gms/lit.</b>			
Casein enzymic hydrolysate		10.00			
Sodium chloride		5.00			
Starch, soluble		1.00			
Phenol red		0.018			
Agar		15.00			
Final pH (at 25°C) : 7.4 ± 0.2					
<b>Directions :</b>					
Suspend 31.02 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add rehydrated contents of 1 vial of Ampicillin Supplement (BF095). Mix well before pouring into sterile Petri plates.					
<b>Principle :</b>					
Casein enzymic hydrolysate in the medium provides essential growth nutrients. Sodium chloride maintains osmotic equilibrium. Ampicillin suppresses the contaminating microflora. Phenol red is the pH indicator.					
<b>QC Tests - (I) Dehydrated Medium</b>					
Colour :		Light yellow to pink			
Appearance :		Homogeneous Free Flowing powder			
<b>(II) Rehydrated medium</b>					
pH (post autoclaving/heating) :		7.4 ± 0.2			
Colour (post autoclaving/heating) :		Red			
Clarity (post autoclaving/heating) :		Clear to slightly opalescent			
<b>(III) Q.C. Test Microbiological</b>					
Cultural characteristics observed after an incubation at 30°C for 24-48 hours with added Ampicillin Supplement (BF095).					
MICROORGANISM (ATCC )		GROWTH	STARCH HYDROLYSIS		
Aeromonas hydrophila (7966)		Luxuriant	Positive, clearing around the colony		
Escherichia coli (25922)		Poor - fair	Negative, no clearing		
Staphylococcus aureus (25923)		Inhibited	Negative, no clearing		
<b>Precautions :</b>	1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
<b>Limitations :</b>	1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
<b>Use :</b>	For isolation, cultivation and differentiation based on starch hydrolysis of Aeromonas hydrophila from foods as per APHA.				
<b>Storage :</b>	Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
<b>Packing :</b>	500 gm. bottle				
<b>Product profile:</b>	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
<b>B704</b>	31g/l	16.129L	7.4 ± 0.2	1 vial of Ampicillin Supplement (BF095)	121°C / 15 minutes

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.