

**TECHNICAL SHEET**

|  |   |                                |           |            |                    |
|--|---|--------------------------------|-----------|------------|--------------------|
| <b>B706</b>  | <b>MILK AGAR WITH CETRIMIDE (TWIN PACK)</b>   |                                |           |            |                    |
| <b>Formula</b>   |   |                                |           |            |                    |
| <b>Ingredients :</b>   |   | <b>gms/lit.</b>                |           |            |                    |
| <b>Part A :</b>  |   |                                |           |            |                    |
| Skim milk Powder   |   | 133.33                         |           |            |                    |
| <b>Part B :</b>  |   |                                |           |            |                    |
| Peptic digest of animal tissue   |   | 3.33                           |           |            |                    |
| Sodium chloride  |   | 1.67                           |           |            |                    |
| Yeast extract  |   | 1.00                           |           |            |                    |
| Cetrimide  |   | 0.40                           |           |            |                    |
| Agar   |   | 20.00                          |           |            |                    |
| Final pH (at 25°C) : 7.3 ± 0.2   |   |                                |           |            |                    |
| <b>Directions :</b>  |   |                                |           |            |                    |
| Suspend 26.4 gms. of Part B in 250 ml. distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 20 minutes. Suspend 133 gms. of Part A in 750 ml of distilled water and sterilize by autoclaving at 15 lbs pressure (121°C) for 5 minutes. After autoclaving mix Part A and B and pour into sterile petri plates.                                     |   |                                |           |            |                    |
| <b>Principle :</b>   |   |                                |           |            |                    |
| Peptic digest of animal tissue, yeast extract and skim milk, provide nitrogen, sulphur, vitamins and other growth nutrients. Sodium chloride maintains osmotic equilibrium. Cetrimide (Cetyl trimethylammonium bromide) is a quaternary ammonium compound which inhibits a wide variety of microorganisms including <i>Pseudomonas</i> species other than <i>Pseudomonas aeruginosa</i> . Agar is the solidifying agent. |   |                                |           |            |                    |
| <b>QC Tests – (I) Dehydrated Medium</b>  |   |                                |           |            |                    |
| Colour :   | Part A - Cream Part B - Light yellow  |                                |           |            |                    |
| Appearance :   | Homogeneous Free Flowing powder   |                                |           |            |                    |
| <b>(II) Rehydrated medium</b>  |   |                                |           |            |                    |
| pH (post autoclaving/heating) :  | 7.3 ± 0.2   |                                |           |            |                    |
| Colour (post autoclaving/heating) :  | Cream   |                                |           |            |                    |
| Clarity (post autoclaving/heating) :   | Opalescent  |                                |           |            |                    |
| <b>(III) Q.C. Test Microbiological</b>   |   |                                |           |            |                    |
| Cultural characteristics observed after 24 - 48 hours at 35 - 37°C.  |   |                                |           |            |                    |
| MICROORGANISM (ATCC)   | GROWTH  | PIGMENT                        |           |            |                    |
| <i>Pseudomonas aeruginosa</i> (27853)  | Good - Luxuriant  | Blue - green                   |           |            |                    |
| <i>Pseudomonas maltophilia</i> (13637)   | Inhibited   | --                             |           |            |                    |
| <i>Escherichia coli</i> (25922)  | Inhibited   | --                             |           |            |                    |
| <b>Precautions :</b>   | 1. For Laboratory Use.<br>2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.          |                                |           |            |                    |
| <b>Limitations :</b>   | 1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium. |                                |           |            |                    |
| <b>Use :</b>   | For detection and enumeration of <i>Pseudomonas aeruginosa</i> in water.  |                                |           |            |                    |
| <b>Storage :</b>   | Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.  |                                |           |            |                    |
| <b>Packing :</b>   | 500 gm. bottle  |                                |           |            |                    |
| <b>Product profile:</b>  | Reconstitution  | Quantity on Preparation (500g) | pH (25°C) | Supplement | Sterilization      |
| <b>B706</b>  | 159.73g/l   | 3.13L                          | 7.3 ± 0.2 | NIL        | 121°C / 15 minutes |

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.