

TECHNICAL SHEET

B752	MAINTENANCE (SCY) MEDIUM				
Formula					
Ingredients :		gms/lit.			
Casein enzymic hydrolysate		0.91			
Yeast extract		0.25			
Sucrose		1.00			
Sodium chloride		0.05			
Dipotassium hydrogen phosphate		0.02			
Potassium iodide		0.03			
Thiamine		0.0004			
Agar		10.00			
Final pH (at 25°C) : 7.3 ± 0.2					
Directions :					
Suspend 12.26 gms. in 1000ml distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15lbs pressure (121°C) for 15 minutes. Aseptically add filter sterilized solution of Cyanocobalamin to a final concentration of 0.01 mg/litre.					
Principle :					
Casein enzymic hydrolysate and yeast extract provide all essential nutrients. Sucrose is the carbohydrate. Sodium chloride maintains osmotic balance. Dipotassium hydrogen phosphate act a buffering agent. Agar is the solidifying agent.					
QC Tests - (I) Dehydrated Medium					
Colour :		Cream			
Appearance :		Homogeneous Free Flowing powder			
(II) Rehydrated medium					
pH (post autoclaving/heating) :		7.3 ± 0.2			
Colour (post autoclaving/heating) :		Cream			
Clarity (post autoclaving/heating) :		Clear to slightly opalescent			
(III) Q.C. Test Microbiological					
Cultural characteristics observed after 48 -72 hrs. at 30°C.					
MICROORGANISM (ATCC)		GROWTH			
Sphaerotilus discophorus		Luxuriant			
Sphaerotilus natans (13338)		Luxuriant			
Thiobacillus thioparas (8158)		Luxuriant			
Precautions :					
1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.					
Limitations :					
1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.					
Use :					
For the maintenance of iron bacteria.					
Storage :					
Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.					
Packing :					
500 gm bottle					
Product profile:					
	Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B752	12.26g/l	40.783L	7.3 ± 0.2	Cyanocobalamin	121°C / 15 minutes

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.