

TECHNICAL SHEET

B773	TERGITOL – 7 BROTH					
Formula						
Ingredients :		gms/lit				
Proteose peptone		5.00				
Yeast extract		3.00				
Lactose		10.00				
Sodium heptadecyl sulphate		0.10				
Bromo thymol blue		0.025				
Final pH (at 25°C) :		6.9 ± 0.2				
Directions :						
Suspend 18.13 gms.in 1000 ml distilled water. Boil to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add 3 ml. of 1% 2, 3, 5, Triphenyl Tetrazolium Chloride (TTC) solution, if desired.						
Principle :						
Proteose peptone provides the carbon and nitrogen sources required for good growth of a wide variety of organisms. Vitamins and cofactors required for growth, as well as additional sources of nitrogen and carbon, are provided by yeast extract.						
QC Tests - (I) Dehydrated Medium						
Colour :		Light to medium yellow				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		6.9 ± 0.2				
Colour (post autoclaving/heating) :		Green				
Clarity (post autoclaving/heating) :		Clear to slightly opalescent				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after 18 – 48 hrs.at 35-37 °C.						
MICROORGANISM (ATCC)		GROWTH		COLOUR OF MEDIUM		
Enterobacter aerogenes (13048)		Luxuriant		yellow		
Escherichia coli (25922)		Luxuriant		yellow		
Salmonella typhimurium (14028)		Luxuriant		Blue		
Shigella flexneri (12022)		Luxuriant		Blue		
Staphylococcus aureus (25923)		Inhibited		-		
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium. 2. Since the medium with TTC permits growth of coliform organisms, this fact must be taken into consideration in the isolation of Candida from specimens. 3. Pour plates do not give satisfactory results. 4. Allow plates to dry with lids slightly ajar for 1-2 hours after dispensing. 5. Reduction of TTC is an irreversible reaction that produces an insoluble formazan compound.				
Use :		For detection and enumeration of coliforms.				
Storage :		Dehydrated medium- below 30°C Prepared medium- Between 2 to 8°C.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B773		18.13 g/l	27.578L	6.9 ± 0.2	1% 2,3,5, Triphenyl Tetrazolium Chloride (TTC) solution	121°C /15 min.

Refer disclaimer Overleaf

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.