

TECHNICAL SHEET

B776		THIOGLYCOLLATE AGAR				
Formula						
Ingredients :		gms/lit.				
Casein enzymic hydrolysate		15.00				
L-Cystine		0.50				
Dextrose		5.50				
Yeast extract		5.00				
Sodium chloride		2.50				
Sodium thioglycollate		0.50				
Resazurin		0.001				
Agar		20.00				
Final pH (at 25°C) :		7.1 ± 0.2				
Directions :						
Suspend 49.0 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Dispense as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.						
Principle :						
Casein enzymic hydrolysate, yeast extract provides nitrogenous compounds, vitamin B and other essential growth nutrients. Dextrose is the fermentable carbohydrate and energy source. Resazurin is the redox indicator. Sodium thioglycollate act as a reducing agent and neutralizes the antibacterial effect of mercurial preservatives and other heavy metal compounds which exert a bacteriostatic effect in the materials under examination. Any increase in the oxygen content is indicated by a colour change of redox indicator, resazurin to red. Agar is the solidifying agent.						
QC Tests - (I) Dehydrated Medium						
Colour :		Cream to yellow				
Appearance :		Homogeneous Free Flowing powder				
(II) Rehydrated medium						
pH (post autoclaving/heating) :		7.1 ± 0.2				
Colour (post autoclaving/heating) :		Light amber				
Clarity (post autoclaving/heating) :		Clear to very slightly opalescent solution				
(III) Q.C. Test Microbiological						
Cultural characteristics observed after 40-48 hrs.at 35-37°C under anaerobic condition						
MICROORGANISM (ATCC)		GROWTH				
Clostridium botulinum (25763)		Luxuriant				
Clostridium perfringenes (12924)		Luxuriant				
Clostridium sporogenes (11437)		Luxuriant				
Precautions :		1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.				
Limitations :		1. Since the nutritional requirements of organisms vary, some strains may be encountered that fail to grow or grow poorly on this medium.				
Use :		For cultivation of anaerobic microorganisms.				
Storage :		Dehydrated medium- below 30°C Prepared medium - At room temperature and away from light.				
Packing :		500 gm. bottle				
Product profile:		Reconstitution	Quantity on Preparation (500g)	pH (25°C)	Supplement	Sterilization
B776		49.0g/l	10.204 L	7.1 ± 0.2	NIL	121°C /15 min.

Refer disclaimer Overleaf

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.