

**BA020**

**NIGROSIN STAIN 10% W/V**

**FORMULA**

**Reagents:**

Nigrosin	10.00 g
Formalin	0.50 ml
Distilled water	100.00 ml

**Directions:**

1. Take a loopful of culture on clear, dry glass slide.
2. Mix it with a loopful of Nigrosin.
3. Mix well and cover with a thin coverslip.
4. Observe under high power or oil immersion lens.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** Negative staining.

**Quality Control:**

**Appearance:** Blackish violet solution

**Microscopy:**

Bacteria : Clear halo surrounding the bacterial cell.  
Background : Dark of carbon particles.

**Storage:** Below 30°C.

**Packing:** 100 ml