

**BA021**

**LACTOPHENOL**

**FORMULA**

**Reagents:**

Phenol crystals	20.00 g
Lactic acid	20.00 ml
Glycerol	20.00 ml
Distilled water	20.00 ml

**Directions :**

1. Take a loopful of fungi culture on clear, dry glass slide.
2. Mix it with a loopful of Lactphenol cotton blue.
3. Tease the fungal culture into thin preparation carefully with the help of wire loop.
4. Place coverslip.
5. Observe under high power.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** Fungal staining

**Quality Control:**

**Appearance:** Pale yellow coloured liquid.

**Microscopy:**

**Storage:** Below 30°C.

**Packing :** 100 ml