

TECHNICAL SHEET

**BA085**

**BARRIT REAGENT B**

**Formula**

**Reagent:**

KOH	40.00 g
Distilled water	100.00 ml

**Appearance:** Colourless solution

**Directions:**

1. Grow test culture in MR-VP (BW543) and incubate at 30-32°C for 24-48 hours.
2. Add 0.2 ml of both reagents (BA084 and BA085) to 10 ml medium after incubation.
3. Shake tubes for 1 minute. Allow to stand for 10-15 minutes.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** Differentiating E. coli from Klebsiella-Enterobacter groups.

**Quality Control:**

**Microbiological testing**

<b>MICROORGANISM (ATCC)</b>	<b>VP TEST</b>
E. coli (25922)	Negative (No red colour formation)
Enterobacter aerogenes (13048)	Positive (Red colour formation)

**Storage:** Below 30°C

**Packing:** 100 ml bottle

**BIOMARK Laboratories-INDIA**

[www.biomarklabs.com](http://www.biomarklabs.com)

**TECHNICAL SHEET**

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 02 of 02