

**BA107 SODIUM HYDROXIDE SOLUTION (1N,0.1N)**

**Formula**

**Preparation:**

**1N NAOH**

Dissolve 40 g of sodium hydroxide in 150 ml of distilled water in a 1000 ml volumetric flask, cool to room temperature, add distilled water slowly to make up final volume to 1000 ml.

**0.1 N NaOH**

Dissolve 4 g of sodium hydroxide in 150 ml of distilled water in a 1000 ml volumetric flask, cool to room temperature, add distilled water slowly to make up final volume to 1000 ml.

**Precautions :**

1. For Laboratory Use.
2. Corrosive.Handle carefully.
3. Exothermic reaction.Cool properly.

**Use:**For volumetric analysis/pH adjustment.

**Quality Control:**

**Appearance:** Colourless liquid

**Storage:** Below 30°C

**Packing :** 500 ml bottle

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.