

BIOMARK Laboratories-INDIA
www.biomarklabs.com
TECHNICAL SHEET

BA122FEHLING SOLUTION NO.2

Formula

Reagent:

Sodium potassium tartrate	176.00 g	
NaOH		77.00 g
Distilled water	500.00 ml	

Appearance: Colourless clear solution

Mix equal volumes of the Fehling solutions No.1 and No.2 at the time of using.

Precautions:

1. For Laboratory Use.

Use: Estimation of inverted sugar.

Storage: Below 30oC

Packing: 500 ml bottle

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.