

**BA128                      STERILE 3.5% POTASSIUM TELLURITE SOLUTION**

**Formula**

**Reagent:**

Potassium tellurite	3.5 g
Distilled water	100.00 ml

**Precautions :**

1. For Laboratory Use.

**Use:** For use in Production of Egg yolk tellurite emulsion.

**Quality Control:**

**Appearance:** Colorless clear liquid.

**Sterility:** Checked for any growth in Nutrient Broth for 24 hrs. at 35-37°C.

**Storage:** Between 2-8°C

**Packing:** 100 ml bottle

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications. The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.