

**BDD011**

**MALTOSE(Ma)**

**SPECIFICATIONS:**

**COLOUR:** White discs of 6 mm diameter with Ma printed on both sides of discs

**DIRECTIONS:**

1. A sugar free medium base Phenol Red Broth Base (B284)/ Andrade Peptone Water (B834) is prepared as per directions given on label.
2. 5 ml of liquid is dispensed with inverted Durham's tube in test tubes and sterilized.
3. A single carbohydrate disc is added to each tube aseptically and inoculated with test organisms.
4. Incubate at 35 - 37°C for 18- 48 hours.

**RESULTS:**

Gas produced during fermentation is collected in the inverted Durham's tube while acid produced changes colour of the medium.

**USE:** It is used for bacterial identification for relevant test organisms for laboratory diagnosis.

**Precautions :**For Laboratory Use.

**Storage :** Store at 2-8°C

**Packing:** 25 discs in vial/cartridge