BIOMARK Laboratories-INDIA www.biomarklabs.com TECHNICAL SHEET

CHLORAMPHENICOL SELECTIVE SUPPLEMENT

BF004

FORMULA	
Ingredients:	
(Per vial sufficient for 500 ml medium)	
Chloramphenicol 50.00 mg	
Directions:	
Wi Ba Mi Ba	ehydrate the contents of 1 vial aseptically with 2 ml of ethanol. Mix ell and aseptically add to 500 ml of sterile, molten SABHI Agar ase(B376) / Rose Bengal Chloramphenicol Agar Base(B307) / Dichloran edium Base with Rose Bengal(B1010) / Aspergillus Differentiation Medium ase(B846)/ DG18 Agar Base(B1439)/DRBC Agar Base (B1438). Mix well nd pour into sterile petri plates.
Precautions :	 For Laboratory Use. Follow proper, established laboratory procedures in handling and disposing of infectious materials.
Use:	An antibiotic supplement recommended for the selective isolation of yeasts and molds.
Storage:	Between 2 to 8°C.
Packing:	5 vials