

BF009

B P SULPHA SUPPLEMENT

FORMULA

Ingredients:

(Per vial sufficient for 1000 ml medium)

Sulphamethazine	50.000 mg
-----------------	-----------

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and add to 1000 ml sterile molten cooled (45-50°C) Baired Parker Agar base (B013) aseptically. Mix well and pour in sterile Petri Plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: This is an antibacterial supplement recommended for suppression of Proteus species growing on Baired Parker Agar base (B013)

Storage: Between 2 to 8°C.

Packing: 5 vials