BF012 BRUCELLA SELECTIVE SUPPLEMENT, MODIFIED

FORMULA

Ingredients:

(per vial sufficient for 500 ml of medium)

Polymixin B	2500 IU
Bacitracin	12500 IU
Nystatin	50000 IU
Cycloheximide	50.00 mg
Nalidixic acid	2.50 mg
Vancomycin	10.00 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 10 ml of 50% methanol. Shake to form uniform suspension. Add the contents to 500 ml of sterile molten media such as Blood agar base no.2(B116) /Blood Agar Base No.2 (w/Agar 1.2%) (B117)/Columbia Blood agar base (B994)/Columbia Blood Agar Base (w/Agar 1%) (B148) with 5-10% v/v inactivated horse serum and 1% w/v sterile dextrose or Brucella Agar base(B125)/ Brucella Broth base(B126) with 5-10% v/v inactivated horse serum. Mix well and pour into sterile Petri plates/tubes.

Precautions : 1. For Laboratory Use.

Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: For selective isolation of Brucella species from milk

Storage: Between 2 to 8°C.

Packing: 1 vial

Disclaimer:

Page 01 of 01

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.