

TECHNICAL SHEET

BF014 CAMPYLOBACTER SUPPLEMENT-II (BUTZLER, MODIFIED)

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Bacitracin	12500.000 IU
Colistin sulphate	5000.000 units
Amphotericin B	5.000 mg
Cephazolin sodium	7.500 mg
Novobiocin	2.500 mg

Directions:

Rehydrate the content of one vial aseptically with 5ml of 50% aqueous Ethanol and mix well to dissolve. Avoid frothing of the solution. Aseptically add the contents to 500 ml of sterile molten cooled (45-50°C) Columbia Blood Agar Base (B994) / Columbia Blood Agar Base w/1% Agar (B148) / Blood Agar Base N0 2 (B116)/Blood Agar Base N0 2 w/1.2% Agar(B117) /Brucella Broth Base (B126)/For 470 ml of Modified Brucella Agar Base (B690) along with 5%v/v defibrinated sheep blood. Mix gently and pour into sterile petri plates.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: Campylobacter supplement recommended for the selective isolation of campylobacter species

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.