

**TECHNICAL SHEET**

**BF017 STREPTO SUPPLEMENT**

**FORMULA**

**Ingredients:**

(Per vial sufficient for 500 ml medium)

Nalidixic Acid	3.75 mg
Neomycin Sulphate	1.06 mg
Polymyxin B Sulphate	8500.00 units

**Directions:** Rehydrate the contents of 1 vial aseptically with 2-5 ml. of sterile distilled water. Aseptically add the contents to 500 ml sterile molten cooled (45-50°C) Columbia Blood Agar Base (B994)/Columbia Blood Agar Base (w/Agar 1%) (B148) /Blood Agar Base No 2 (B116)/Blood Agar Base No.2 (w/Agar 1.2%) (B117). The medium may be enriched with 5-7% v/v Sterile defibrinated blood. Mix gently and pour into sterile Plates

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** Strepto Supplement recommended for the selective isolation of Streptococcus species

**Storage:** Between 2 to 8°C.

**Packing:** 5 vials

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.