

TECHNICAL SHEET

BF020	SULPHA SUPPLEMENT
FORMULA	
Ingredients:	
(per vial sufficient for 500 ml medium)	
Sodium sulpha acetamide	500 mg
Sodium Mandelate	125 mg
Directions:	
Rehydrated the content of 1 vial with 5.0 ml sterile distilled water. Aseptically add these to 500 ml sterile, molten, cooled (45-50°C) Brilliant green agar, Modified(B123)/ Brilliant green agar (w/Agar 1.2%) (B124)/Brilliant green agar w/phosphates (B439, BW439)/ Brilliant green agar w/phosphates (B1504). Mix well and pour in sterile petri plates.	
Precautions :	1. For Laboratory Use. 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.
Use:	An antibacterial supplement recommended for selective isolation of Salmonella species.
Storage:	Between 2 to 8°C.
Packing:	5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.