BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF021

HORSE SERUM

FORMULA

Ingredient:

Horse serum

Appearance: Clear liquid, may have slight particles

Colour: Straw to amber

Directions:

For Brucella agar base(B125), Brucella broth base(B126), Fletcher Leptospira Medium Base(B1039), Frey Mycoplasma broth base(B529), Mycoplasma agar base (B265) Mycoplasma broth base w/CV (B266), Mycoplasma broth base w/CV (B267), Mycoplasma Synoviae Medium(B555), Trichomonas Agar Base(B337), Trichomonas Broth BaseNo.2(B075), Mycoplasma Urogenital Broth Base (B1233), GBS Medium Base (B1048), Islam's Medium Base for Gr. B Streptococci (B1073), Krischner Medium Base, Modified (B1087)

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.
- 3. Thawing of the serum should be done as quickly as possible.

Use: For selective isolation & cultivation of Brucella, Campylobacter, Leptospira Mycoplasma and Trichomonas species.

Storage: Store at -20°C

Packing: 100 ml

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01