### **BIOMARK Laboratories-INDIA**

### www.biomarklabs.com

# **TECHNICAL SHEET**

### BF026

# LEGIONELLA SUPPLEMENT (TWIN PACK)

### FORMULA Ingredients:

(Per vial, sufficient for 500 ml medium)

# Part A

I UICA	
L-Cysteine hydrochloride	200 mg
Part B	
Ferric pyrophosphate, soluble	125mg
Distilled water	5ml

# **Directions:**

Rehydrate the contents of part A with 5 ml sterile distilled water. Warm up the refrigerated solution of Part B to room temperature and aseptically add Part A and Part B to 490 ml of sterile, molten, cooled (45-50°C) Legionella Agar Base (B221) / Buffered Charcoal Yeast Extract Agar Base (B127) / Modified Buffered Charcoal Agar Base (B627) along with Legionella Selective Supplement as recommended. Mix well and pour into sterile Petri plates

# **Precautions :**

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** Recommended for enhanced growth of Legionella species

Storage: Between 2 to 8°C.

Packing :5 vial

#### Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Page 01 of 01