BIOMARK Laboratories-INDIA

www.biomarklabs.com

TECHNICAL SHEET

BF031 CPC SUPPLEMENT

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Colistin methane sulphonate
Polymyxin B sulphate

6,80,000 Units 50,000 Units

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml. sterile distilled water. Mix well and aseptically add to 500 ml. of sterile, molten CPC Agar Base(B960). Mix well and pour into sterile petri plates.

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antibiotic supplement recommended for the selective isolation and numeration of Vibrio species from foods.

Storage: Between 2 to 8°C.

Packing: 5 vial