

**TECHNICAL SHEET**

**BF032**

**PARK AND SANDERS SELECTIVE SUPPLEMENT II**

**FORMULA**

**Ingredients:**

(Per vial)

Cefoperazone	32.00 mg
Amphotericin B	5.00 mg

**Directions:**

Rehydrate the contents of 1 vial aseptically with 5 ml. of 50% ethanol Mix well and aseptically add to 995 ml of sterile, Park and Sanders Broth Base(B1480)/Campylobacter Cefex agar base(B965) containing 50ml sterile, defibrinated blood and Park and Sanders Selective Supplement (BF087). Mix well and dispense.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** An antibiotic supplement recommended for selective isolation and enumeration of thermo tolerant Campylobacter species

**Storage:** Between 2 to 8°C.

**Packing:** 5 vials

**Disclaimer:**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.