#### **BIOMARK Laboratories-INDIA**

#### www.biomarklabs.com

#### **TECHNICAL SHEET**

## BF033 CAMPYLOBACTER SELECTIVE SUPPLEMENT IV (PRESTON)

## **FORMULA**

# **Ingredients:**

(Per vial, sufficient for 500 ml medium)

Polymyxin B sulphate 2500 IU
Rifampicin 5 mg
Trimethoprim 5 mg
Cycloheximide 50 mg

## **Directions:**

Rehydrate the contents of 1 vial aseptically with 5 ml. of 50% Acetone. Mix well and aseptically add it to 500 ml sterile molten Campylobacter Enrichment Broth Base (B966)/Preston Agar Base(B1182) along with 25 ml sterile, lysed horse blood. Mix well and dispense into sterile tubes or Petri plates.

**Precautions:** 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and

disposing of infectious materials.

**Use:** This supplement recommended for the selective isolation of

Campylobacter species.

**Storage:** Between 2 to 8°C.

**Packing:** 5 vials