

BF033 CAMPYLOBACTER SELECTIVE SUPPLEMENT IV (PRESTON)

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Polymyxin B sulphate 2500 IU
Rifampicin 5 mg
Trimethoprim 5 mg
Cycloheximide 50 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml. of 50% Acetone. Mix well and aseptically add it to 500 ml sterile molten Campylobacter Enrichment Broth Base (B966)/Preston Agar Base(B1182) along with 25 ml sterile, lysed horse blood. Mix well and dispense into sterile tubes or Petri plates.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: This supplement recommended for the selective isolation of Campylobacter species.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARK LABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.