

BF036

HAEMOGLOBIN POWDER

FORMULA

Ingredients:

(Per bottle)

Haemoglobin Powder 50 g/100 g

Directions:

A specially prepared powder whose 2% w/v aqueous solution is autoclavable. The aqueous solution is chocolate brown, opaque and contains flocculent dispersible precipitate. It is used for 500 ml medium preparation of GC Agar Base (B1049) - 5 g / Thayer Martin Medium Base (B333) - 5 g / Tellurite Blood Agar Base (B085) - 10 g/ Chocolate Agar Base (B980) -5 g / Modified Protease Agar - 5 g / Transgrow Medium Base (B078) - 2 g.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An enrichment supplement recommended for isolation of Neisseria.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.