BIOMARK Laboratories-INDIA www.biomarklabs.com

TECHNICAL SHEET

BF040 G.VAGINALIS SELECTIVE SUPPLEMENT

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Gentamicin sulphate 2 mg Nalidixic acid 15 mg Amphotericin B 1 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Columbia Blood Agar Base (B148) along with 5% v/v sterile defibrinated sheep blood. Mix well and pour into sterile petri plates.

Precautions:

- 1. For Laboratory Use.
- 2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use :An antibiotic supplement recommended for the selective isolation of Gardnerella vaginalis .

Storage: Between 2 to 8°C.

Packing: 5 vial